



3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone 800-325-5832 • (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sia.com  
sigma-aldrich.com

## Product Information

### Cytochrome c from bovine heart

Product Number **C 3131**  
Storage Temperature -0 °C

#### Product Description

CAS Number: 9007-43-6  
Molecular Weight: 12.23 kDa<sup>1</sup>  
pI: 10.37-10.80<sup>2</sup>  
 $\lambda_{\text{max}}$ : 550 nm (reduced)  
Extinction coefficient:  $E^{\text{mM}} = 28.0$  (reduced).

Cytochrome c is an electron-carrying mitochondrial protein. It is a small heme protein containing a single polypeptide chain and a single heme group, which is covalently attached to the polypeptide. The ready fluctuation of cytochrome c within the cell between ferrous and ferric states, makes it an efficient biological electron-transporter and it plays a vital role in cellular oxidations in both plants and animals. It is generally regarded as a universal catalyst of respiration, forming the essential electron-bridge between the respirable substrates and oxygen.

This product is mainly the oxidized form of cytochrome c, often referred to as ferricytochrome c. The preparation of reduced cytochrome c can be made utilizing either sodium dithionite or sodium ascorbate followed by gel filtration.<sup>3,4</sup>

#### Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

#### Preparation Instructions

Cytochrome c is soluble in water (10 mg/ml), yielding a dark red to red brown solution.

Cytochrome c solutions can be prepared in 50 mM phosphate buffer at neutral pH. Solutions can be stored frozen as single-use aliquots at -20 °C for long term storage or at 2-8 °C for approximately 1-2 weeks.

#### References

1. Nakashima, T., et al., The amino acid sequence of bovine heart cytochrome c. *J. Biol. Chem.*, **241(5)**, 1166-1177 (1966).
2. Righetti, P.G., and Caravaggio, T., J. Chromatography, **127**, 1-28 (1976).
3. Dixon, H.B., and McIntosh, R., Reduction of methaemoglobin in haemoglobin samples using gel filtration for continuous removal of reaction products. *Nature*, **213(74)**, 399-400 (1967).
4. Errede, B., et al., Preparation and properties of complex IV (ferrocytochrome c:Oxygen Oxidoreductase EC 1.9.3.1). *Methods in Enzymology*, **53D**, 40-47 (1978).

ARO/NSB 8/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.